

SIX3 Polyclonal Antibody

Cat No: HR1AP10851

For research use only

Overview

Product Name	SIX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	SIX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SIX homeobox 3(SIX3) Homo sapiens This gene encodes a member of the sine oculis homeobox transcription factor family. The encoded protein plays a role in eye development. Mutations in this gene have been associated with holoprosencephaly type 2. [provided by RefSeq, Oct 2009],