

NKX22 Polyclonal Antibody

Cat No: HR1AP10847

For research use only

Overview

Product Name	NKX22 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	NKX2-2 NKX2.2 NKX2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NK2 homeobox 2(NKX2-2) Homo sapiens The protein encoded by this gene contains a homeobox domain and may be involved in the morphogenesis of the central nervous system. This gene is found on chromosome 20 near NKX2-4, and these two genes appear to be duplicated on chromosome 14 in the form of TITF1 and NKX2-8. The encoded protein is likely to be a nuclear transcription factor. [provided by RefSeq, Jul 2008],