

NOBOX Polyclonal Antibody

Cat No: HR1AP10850

For research use only

Overview

Product Name	NOBOX Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	NOBOX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NOBOX oogenesis homeobox(NOBOX) Homo sapiens This homeobox gene encodes a transcription factor that is thought to play a role in oogenesis. In mice, it is essential for folliculogenesis and regulation of oocyte-specific genes. Defects in this gene result in premature ovarian failure type 5.[provided by RefSeq, May 2011],