

NKX23 Polyclonal Antibody

Cat No: HR1AP10848

For research use only

Overview

Product Name	NKX23 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	NKX2-3 NKX23 NKX2C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NK2 homeobox 3(NKX2-3) Homo sapiens This gene encodes a homeodomain-containing transcription factor. The encoded protein is a member of the NKX family of homeodomain transcription factors. Studies of similar proteins in mouse and rat have indicated a potential role in cellular differentiation.[provided by RefSeq, Mar 2010],