

# MEIS3 Polyclonal Antibody

Cat No: HR1AP10846

For research use only

## Overview

Product Name	MEIS3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	MEIS3 MRG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Meis homeobox 3(MEIS3) Homo sapiens This gene encodes a homeobox protein and probable transcriptional regulator. The orthologous protein in mouse controls expression of 3-phosphoinositide dependent protein kinase 1, which promotes survival of pancreatic beta-cells. [provided by RefSeq, Sep 2016],