

DBX1 Polyclonal Antibody

Cat No: HR1AP10827

For research use only

Overview

Product Name	DBX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	DBX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>DBX1 (Developing Brain Homeobox 1) is a Protein Coding gene. Diseases associated with DBX1 include Central Hypoventilation Syndrome, Congenital. Could have a role in patterning the central nervous system during embryogenesis. Has a key role in regulating the distinct phenotypic features that distinguish two major classes of ventral interneurons, V0 and V1 neurons. Regulates the transcription factor profile, neurotransmitter phenotype, intraspinal migratory path and axonal trajectory of V0 neurons, features that differentiate them from an adjacent set of V1 neurons (By similarity).</p>