

CUX2 Polyclonal Antibody

Cat No: HR1AP10826

For research use only

Overview

Product Name	CUX2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	163kDa
GeneID?Human?	CUX2 CUTL2 KIAA0293
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cut like homeobox 2(CUX2) Homo sapiens This gene encodes a protein which contains three CUT domains and a homeodomain; both domains are DNA-binding motifs. A similar gene, whose gene product possesses different DNA-binding activities, is located on chromosome on chromosome 7. Two pseudogenes of this gene have been identified on chromosomes 10 and 4. [provided by RefSeq, Jan 2013],