

HESX1 Polyclonal Antibody

Cat No: HR1AP10825

For research use only

Overview

Product Name	HESX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	HESX1 HANF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	HESX homeobox 1(HESX1) Homo sapiens This gene encodes a conserved homeobox protein that is a transcriptional repressor in the developing forebrain and pituitary gland. Mutations in this gene are associated with septo-optic dysplasia, HESX1-related growth hormone deficiency, and combined pituitary hormone deficiency. [provided by RefSeq, Jul 2008],