

HBXIP Polyclonal Antibody

Cat No: HR1AP10818

For research use only

Overview

Product Name	HBXIP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	HBXIP XIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	late endosomal/lysosomal adaptor, MAPK and MTOR activator 5(LAMTOR5) Homo sapiens This gene encodes a protein that specifically complexes with the C-terminus of hepatitis B virus X protein (HBx). The function of this protein is to negatively regulate HBx activity and thus to alter the replication life cycle of the virus. [provided by RefSeq, Jul 2008],