

# HEAT1 Polyclonal Antibody

Cat No: HR1AP10817

For research use only

## Overview

Product Name	HEAT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	235kDa
GeneID?Human?	HEATR1 BAP28
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in nucleolar processing of pre-18S ribosomal RNA. Involved in ribosome biosynthesis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the HEATR1/UTP10 family.,similarity:Contains 1 HEAT repeat.,