

GNL3 Polyclonal Antibody

Cat No: HR1AP10815

For research use only

Overview

Product Name	GNL3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	GNL3 E2IG3 NS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein nucleolar 3(GNL3) Homo sapiens The protein encoded by this gene may interact with p53 and may be involved in tumorigenesis. The encoded protein also appears to be important for stem cell proliferation. This protein is found in both the nucleus and nucleolus. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],