

## MSS4 Polyclonal Antibody

Cat No: HR1AP10808

For research use only

### Overview

Product Name	MSS4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	RABIF MSS4 RASGRF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RAB interacting factor(RABIF) Homo sapiens This gene encodes a member of the SCE4/YPT1/RAB family of small GTP-binding proteins that are involved in the regulation of intracellular vesicular transport. This protein stimulates GTP-GDP exchange in SEC4, and to a lesser extent in YPT1 and RAB3A, and may play a general role in vesicular transport. [provided by RefSeq, Oct 2011],