

MCF2L Polyclonal Antibody

Cat No: HR1AP10807

For research use only

Overview

Product Name	MCF2L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	MCF2L KIAA0362
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MCF.2 cell line derived transforming sequence like(MCF2L) Homo sapiens This gene encodes a guanine nucleotide exchange factor that interacts specifically with the GTP-bound Rac1 and plays a role in the Rho/Rac signaling pathways. A variant in this gene was associated with osteoarthritis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],