

# REM2 Polyclonal Antibody

Cat No: HR1AP10804

For research use only

## Overview

Product Name	REM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	REM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:It is uncertain whether Met-1 or Met-59 is the initiator.,function:Binds GTP saturably and exhibits a low intrinsic rate of GTP hydrolysis.,similarity:Belongs to the small GTPase superfamily. RGK family.,