

# RIT2 Polyclonal Antibody

Cat No: HR1AP10805

For research use only

## Overview

Product Name	RIT2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RIT2 RIN ROC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Ras like without CAAX 2(RIT2) Homo sapiens RIN belongs to the RAS (HRAS; MIM 190020) superfamily of small GTPases (Shao et al., 1999 [PubMed 10545207]).[supplied by OMIM, Mar 2008],