

GEM Polyclonal Antibody

Cat No: HR1AP10802

For research use only

Overview

Product Name	GEM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	GEM KIR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GTP binding protein overexpressed in skeletal muscle(GEM) Homo sapiens The protein encoded by this gene belongs to the RAD/GEM family of GTP-binding proteins. It is associated with the inner face of the plasma membrane and could play a role as a regulatory protein in receptor-mediated signal transduction. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008].