

# REM1 Polyclonal Antibody

Cat No: HR1AP10803

For research use only

## Overview

Product Name	REM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	REM1 GES REM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RRAD and GEM like GTPase 1(REM1) Homo sapiens The protein encoded by this gene is a GTPase and member of the RAS-like GTP-binding protein family. The encoded protein is expressed in endothelial cells, where it promotes reorganization of the actin cytoskeleton and morphological changes in the cells. [provided by RefSeq, Jul 2008],