

## DIRA3 Polyclonal Antibody

Cat No: HR1AP10801

For research use only

### Overview

Product Name	DIRA3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	DIRAS3 ARHI NOEY2 RHOI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>DIRAS family GTPase 3(DIRAS3) Homo sapiens This gene encodes a member of the ras superfamily. This gene is imprinted gene with monoallelic expression of the paternal allele which is associated with growth suppression. The encoded protein acts as a tumor suppressor whose function is abrogated in many ovarian and breast cancers. This protein may also play a role autophagy in certain cancer cells by regulating the autophagosome initiation complex. [provided by RefSeq, Nov 2015],</p>