



DIRA2 Polyclonal Antibody

Cat No: HR1AP10800

For research use only

Overview

Product Name	DIRA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	DIRAS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DIRAS family GTPase 2(DIRAS2) Homo sapiens DIRAS2 belongs to a distinct branch of the functionally diverse Ras (see HRAS; MIM 190020) superfamily of monomeric GTPases.[supplied by OMIM, Apr 2004],