

GIMA7 Polyclonal Antibody

Cat No: HR1AP10796

For research use only

Overview

Product Name	GIMA7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	GIMAP7 IAN7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GTPase, IMAP family member 7(GIMAP7) Homo sapiens This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008],