

## GTPB1 Polyclonal Antibody

Cat No: HR1AP10798

For research use only

### Overview

Product Name	GTPB1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	GTPBP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GTP binding protein 1(GTPBP1) Homo sapiens This gene is upregulated by interferon-gamma and encodes a protein that is a member of the AGP11/GTPBP1 family of GTP-binding proteins. A structurally similar protein has been found in mouse, where disruption of the gene for that protein had no observable phenotype. [provided by RefSeq, Jul 2008],