

RGPA1 Polyclonal Antibody

Cat No: HR1AP10797

For research use only

Overview

Product Name	RGPA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	223kDa
GeneID?Human?	RALGAPA1 GARNL1 KIAA0884 TULIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Ral GTPase activating protein catalytic alpha subunit 1(RALGAPA1) Homo sapiens This gene encodes a major subunit of the RAL-GTPase activating protein. A similar protein in mouse binds E12, a transcriptional regulator of immunoglobulin genes. The mouse protein also functions in skeletal muscle by binding to the regulatory 14-3-3 proteins upon stimulation with insulin or muscle contraction. A pseudogene of this gene has been identified on chromosome 9. [provided by RefSeq, Oct 2016],</p>