

## GASP2 Polyclonal Antibody

Cat No: HR1AP10791

For research use only

### Overview

Product Name	GASP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	GPRASP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor associated sorting protein 2(GPRASP2) Homo sapiens The protein encoded by this gene is a member of a family that regulates the activity of G protein-coupled receptors (GPCRs). The encoded protein has been shown to be capable of interacting with several GPCRs, including the M1 muscarinic acetylcholine receptor and the calcitonin receptor. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010],