

GON4L Polyclonal Antibody

Cat No: HR1AP10790

For research use only

Overview

Product Name	GON4L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	246kDa
GeneID?Human?	GON4L GON4 KIAA1606
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Enhances transcription activation by YY1.,sequence caution:Translated as Tyr.,similarity:Contains 1 Myb-like domain.,similarity:Contains 2 PAH (paired amphipathic helix) repeats.,subunit:Interacts with YY1.,tissue specificity:Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels in healthy liver. Highly expressed in most liver carcinomas.,</p>