

GDN Polyclonal Antibody

Cat No: HR1AP10785

For research use only

Overview

Product Name	GDN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	SERPINE2 P17 PN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>serpin family E member 2(SERPINE2) Homo sapiens This gene encodes a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. Thrombin, urokinase, plasmin and trypsin are among the proteases that this family member can inhibit. This gene is a susceptibility gene for chronic obstructive pulmonary disease and for emphysema. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],</p>