

CXB5 Polyclonal Antibody

Cat No: HR1AP10781

For research use only

Overview

Product Name	CXB5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	GJB5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	gap junction protein beta 5(GJB5) Homo sapiens This gene encodes a member of the beta-type (group I) connexin family. The encoded protein is a gap junction protein involved in intercellular communication related to epidermal differentiation and environmental sensing. This gene has been linked to non-small cell lung cancer. [provided by RefSeq, Nov 2012],