

GCFC1 Polyclonal Antibody

Cat No: HR1AP10782

For research use only

Overview

Product Name	GCFC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	GCFC1 C21orf66 GCFC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PAX3 and PAX7 binding protein 1(PAXBP1) Homo sapiens This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009],