

## GCP2 Polyclonal Antibody

Cat No: HR1AP10780

For research use only

### Overview

Product Name	GCP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GeneID?Human?	TUBGCP2 GCP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>alternative products:Additional isoforms may exist. May be produced at very low levels due to a premature stop codon in the mRNA, leading to nonsense-mediated mRNA decay,function:Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.,similarity:Belongs to the GCP family.,subunit:Gamma-tubulin complex is composed of gamma-tubulin, GCP2, GCP3, GCP4, GCP5 and GCP6.,tissue specificity:Ubiquitously expressed.,</p>