

## GBRR1 Polyclonal Antibody

Cat No: HR1AP10778

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | GBRR1 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,ELISA  |
| Species Reactivity    | Human   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 52kDa   |
| GeneID?Human?         | GABRR1  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     |   |
| Background            | <p>gamma-aminobutyric acid type A receptor rho1 subunit(GABRR1) Homo sapiens GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. GABRR1 is a member of the rho subunit family. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012],</p> |