

FX4L5 Polyclonal Antibody

Cat No: HR1AP10774

For research use only

Overview

Product Name	FX4L5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	FOXD4L5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	miscellaneous:5 FOXD4-like proteins (FOXD4L2, FOXD4L3, FOXD4L4, FOXD4L5 and FOXD4L6) are encoded by a strongly repeated region of the 9q12 chromosome region. They are very similar and it is therefore difficult to provide a clear and unambiguous protein sequence. Our sequences are in agreement with HGNC nomenclature.,similarity:Contains 1 fork-head DNA-binding domain.,