

FIGL1 Polyclonal Antibody

Cat No: HR1AP10771

For research use only

Overview

Product Name	FIGL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	FIGL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:ATP + H(2)O = ADP + phosphate.,cofactor:Magnesium.,function:May regulate osteoblast proliferation and differentiation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the AAA ATPase family.,subunit:Hexamer.,