

ETAA1 Polyclonal Antibody

Cat No: HR1AP10768

For research use only

Overview

Product Name	ETAA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	ETAA1 ETAA16
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,tissue specificity:Expressed at high levels in the brain, liver kidney and Ewing tumor cell lines.,