

DHB12 Polyclonal Antibody

Cat No: HR1AP10763

For research use only

Overview

Product Name	DHB12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	HSD17B12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hydroxysteroid 17-beta dehydrogenase 12(HSD17B12) Homo sapiens This gene encodes a very important 17beta-hydroxysteroid dehydrogenase (17beta-HSD) that converts estrone into estradiol in ovarian tissue. This enzyme is also involved in fatty acid elongation. [provided by RefSeq, Oct 2011],