

EDEM3 Polyclonal Antibody

Cat No: HR1AP10760

For research use only

Overview

Product Name	EDEM3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	EDEM3 C1orf22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ER degradation enhancing alpha-mannosidase like protein 3(EDEM3) Homo sapiens Quality control in the endoplasmic reticulum (ER) ensures that only properly folded proteins are retained in the cell through recognition and degradation of misfolded or unassembled proteins. EDEM3 belongs to a group of proteins that accelerate degradation of misfolded glycoproteins in the ER (Hirao et al., 2006 [PubMed 16431915]).[supplied by OMIM, Mar 2008].