

NTHL1 Polyclonal Antibody

Cat No: HR1AP10756

For research use only

Overview

Product Name	NTHL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	NTHL1 NTH1 OCTS3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nth like DNA glycosylase 1(NTHL1) Homo sapiens The protein encoded by this gene is a DNA N-glycosylase of the endonuclease III family. Like a similar protein in E. coli, the encoded protein has DNA glycosylase activity on DNA substrates containing oxidized pyrimidine residues and has apurinic/aprymidinic lyase activity. [provided by RefSeq, Oct 2008],