

## ENDD1 Polyclonal Antibody

Cat No: HR1AP10755

For research use only

### Overview

Product Name	ENDD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	ENDOD1 KIAA0830
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May act as a DNase and a RNase .,similarity:Belongs to the DNA/RNA non-specific endonuclease family.,