

ELP3 Polyclonal Antibody

Cat No: HR1AP10754

For research use only

Overview

Product Name	ELP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	ELP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	elongator acetyltransferase complex subunit 3(ELP3) Homo sapiens ELP3 is the catalytic subunit of the histone acetyltransferase elongator complex, which contributes to transcript elongation and also regulates the maturation of projection neurons (Creppe et al., 2009 [PubMed 19185337]).[supplied by OMIM, Apr 2009],