

EF1A2 Polyclonal Antibody

Cat No: HR1AP10751

For research use only

Overview

Product Name	EF1A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	EEF1A2 EEF1AL STN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	eukaryotic translation elongation factor 1 alpha 2(EEF1A2) Homo sapiens This gene encodes an isoform of the alpha subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This isoform (alpha 2) is expressed in brain, heart and skeletal muscle, and the other isoform (alpha 1) is expressed in brain, placenta, lung, liver, kidney, and pancreas. This gene may be critical in the development of ovarian cancer. [provided by RefSeq, Mar 2014],