

UBR2 Polyclonal Antibody

Cat No: HR1AP10749

For research use only

Overview

Product Name	UBR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	193kDa
GeneID?Human?	UBR2 C6orf133 KIAA0349
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin protein ligase E3 component n-recognin 2(UBR2) Homo sapiens This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.[provided by RefSeq, May 2010],