

RN123 Polyclonal Antibody

Cat No: HR1AP10744

For research use only

Overview

Product Name	RN123 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	144kDa
GeneID?Human?	RNF123 KPC1 FP1477
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ring finger protein 123(RNF123) Homo sapiens The protein encoded by this gene contains a C-terminal RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions, and an N-terminal SPRY domain. This protein displays E3 ubiquitin ligase activity toward the cyclin-dependent kinase inhibitor 1B which is also known as p27 or KIP1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],