

RNF41 Polyclonal Antibody

Cat No: HR1AP10742

For research use only

Overview

Product Name	RNF41 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	RNF41 FLRF NRDP1 SBB103
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ring finger protein 41(RNF41) Homo sapiens This gene encodes an E3 ubiquitin ligase. The encoded protein plays a role in type 1 cytokine receptor signaling by controlling the balance between JAK2-associated cytokine receptor degradation and ectodomain shedding. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2011],