

MARH6 Polyclonal Antibody

Cat No: HR1AP10738

For research use only

Overview

Product Name	MARH6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	MARCH6 KIAA0597 RNF176 TEB4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>membrane associated ring-CH-type finger 6(MARCH6) Homo sapiens This gene encodes a member of a family of membrane-associated E3 ubiquitin ligases containing RING-CH-type zinc finger motifs. Ubiquitination of type II deiodinase by the encoded protein is an important regulatory step in thyroid hormone signalling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],</p>