

LRSAM1 Polyclonal Antibody

Cat No: HR1AP10737

For research use only

Overview

Product Name	LRSAM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	LRSAM1 TAL UNQ6496/PRO21356
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>leucine rich repeat and sterile alpha motif containing 1(LRSAM1) Homo sapiens This gene encodes a ring finger protein involved in a variety of functions, including regulation of signaling pathways and cell adhesion, mediation of self-ubiquitylation, and involvement in cargo sorting during receptor endocytosis. Mutations in this gene have been associated with Charcot-Marie-Tooth disease. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jan 2012],</p>