

HACE1 Polyclonal Antibody

Cat No: HR1AP10734

For research use only

Overview

Product Name	HACE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GeneID?Human?	HACE1 KIAA1320
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>HECT domain and ankyrin repeat containing E3 ubiquitin protein ligase 1(HACE1) Homo sapiens This gene encodes a HECT domain and ankyrin repeat-containing ubiquitin ligase. The encoded protein is involved in specific tagging of target proteins, leading to their subcellular localization or proteasomal degradation. The protein is a potential tumor suppressor and is involved in the pathophysiology of several tumors, including Wilms tumor. [provided by RefSeq, Mar 2016].</p>