

HECD1 Polyclonal Antibody

Cat No: HR1AP10735

For research use only

Overview

Product Name	HECD1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	287kDa
GeneID?Human?	HECTD1 KIAA1131
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the UPL family. K-HECT subfamily.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,similarity:Contains 1 MIB/HERC2 domain.,similarity:Contains 4 ANK repeats.,