

## DYLT3 Polyclonal Antibody

Cat No: HR1AP10730

For research use only

### Overview

Product Name	DYLT3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	DYNLT3 TCTE1L TCTE1XL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dynein light chain Tctex-type 3(DYNLT3) Homo sapiens This gene encodes a member of a subclass of dynein light chains. The encoded protein homodimerizes and forms the light chain component of the cytoplasmic dynein motor protein complex. This protein may be important for binding dynein to specific cargos including the spindle checkpoint protein BUB3. This protein may also function independently of dynein as a transcriptional modulator. Pseudogenes of this gene are found on chromosomes 2 and 20.[provided by RefSeq, Mar 2010],</p>