

DLRB1 Polyclonal Antibody

Cat No: HR1AP10728

For research use only

Overview

Product Name	DLRB1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	DYNLRB1 BITH DNCL2A DNLC2A ROBLD1 HSPC162
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dynein light chain roadblock-type 1(DYNLRB1) Homo sapiens This gene is a member of the roadblock dynein light chain family. The encoded cytoplasmic protein is capable of binding intermediate chain proteins, interacts with transforming growth factor-beta, and has been implicated in the regulation of actin modulating proteins. Upregulation of this gene has been associated with hepatocellular carcinomas, suggesting that this gene may be involved in tumor progression. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 12 and 18. [provided by RefSeq, Aug 2013].</p>