

RPAB3 Polyclonal Antibody

Cat No: HR1AP10725

For research use only

Overview

Product Name	RPAB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	POLR2H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RNA polymerase II subunit H(POLR2H) Homo sapiens The three eukaryotic RNA polymerases are complex multisubunit enzymes that play a central role in the transcription of nuclear genes. This gene encodes an essential and highly conserved subunit of RNA polymerase II that is shared by the other two eukaryotic DNA-directed RNA polymerases, I and III. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013],</p>