

RPAC2 Polyclonal Antibody

Cat No: HR1AP10723

For research use only

Overview

Product Name	RPAC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	POLR1D
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RNA polymerase I subunit D(POLR1D) Homo sapiens The protein encoded by this gene is a component of the RNA polymerase I and RNA polymerase III complexes, which function in the synthesis of ribosomal RNA precursors and small RNAs, respectively. Mutations in this gene are a cause of Treacher Collins syndrome (TCS), a craniofacial development disorder. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2011],</p>