



RPC9 Polyclonal Antibody

Cat No: HR1AP10721

For research use only

Overview

Product Name	RPC9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	CRCP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CGRP receptor component(CRCP) Homo sapiens This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],